EastGene Laboratories Inc
42-18 162nd Street, Flushing, NY 11358
Phone: (718) 463 4888, Fax: (718) 463 4889
Email: service@eastgene.com

Saliva Collection Instruction

INTEND USE:

This instruction is intended for the collection of the normal saliva that naturally pools into the mouth. No coughing or sniffing prior to sample collection is required.

STEP BY STEP COLLECTION:

1. Water should be avoided 10 minutes prior to collection. Other drinks, food, and nasal sprays should be avoided for half an hour before sample collection.
2. Before collection, wash your hands thoroughly for 20 seconds and dry.
3. Before collection, ensure collection tube is labelled with the correct identifying information, including patient name, date of birth and collection date.
4. Remove the lid of the wide-mouth sterile collection container (for example: 50 mL tube or urine cup), and begin pooling saliva in your mouth, and then gently expel your saliva into the container until at least 1.0 mL has been collected (see the image with black shading as below). The total volume measured is to exclude any bubbles.

Tip: Saliva production can be stimulated by thinking about food (favorite foods, upcoming meals, etc.) or about the saliva collection itself.

5. Once at least 1.0 mL has been collected, securely replace the lid of the collection container.
6. Sterilize the container surface with 70% ethanol or a disinfecting wipe, and place the sample in an appropriately labeled biohazard bag.

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